

Certificate of Authenticity

Full Cast Crown

The dental restorations used for your case were constructed of FDA approved 510k materials and hand crafted in an ISO 13485:2003 certified dental laboratory.
This is your assurance of quality and safety in materials and construction.

Item	Product Name	Manufacturer	Country of Origin	FDA 510k	Elements & Compounds
High Noble Alloy	Academy	Ivoclar	Liechtenstein	K002872	Au 77%, Ag 12.7%, Cu 8.5%, Pt, In, Ir, Zn, Ta <1%
High Noble Alloy	XL-X	Ivoclar	Liechtenstein	K061851	Au 62.8%, Ag 16.1%, Cu 15.7%, Pd 3.9%, Ir & Zn <1%
Noble Alloy	Minigold	Ivoclar	Liechtenstein	K905326	Au 40.0 %, Ag 47.0 %, Cu 7.5 %, Pd 4.0 %, Zn 1.0 %, < 1.0 % In, Ir
Non-Precious Alloy	Wirobond C	Bego	Germany	K032136	Co 61%, Cr 26%, Mo 6%, W 5%, Si, Fe, Ce each < 2%

**** This product is made entirely from one of the above alloys ****



MODERN
Dental Laboratory USA

Heikas Martirosian

Heikas Martirosian, CDT
Technical Director

