

Certificate of Authenticity

Shaded Cercon® Zirconia

The dental restorations used for your case were constructed of FDA approved 510k materials and hand crafted in an ISO 13485:2003 certified dental laboratory.

This is your assurance of quality and safety in materials and construction.

| Item | Product Name | Manufacturer | Country of Origin | FDA 510k | Elements & Compounds |
|----------------------|----------------|--------------|-------------------|----------|--|
| Zirconium oxide core | Cercon | Dentsply | Germany | K013230 | ZrO ₂ 93%, Y ₂ O ₃ 5%, Hf ₂ O ₃ <2%, other oxides <1% |
| Coloring Liquid | ZirkonZahn Ice | ZirkonZahn | Germany | K061851 | Colored Salts: PEG, FeCl ₃ , MnCl ₂ , CrCl ₂ |
| Porcelain | Vintage ZR | Shofu Dental | USA | K052033 | Ceramic powder made from SiO ₂ , Al ₂ O ₃ , K ₂ O, Na ₂ O, CaO, B ₂ O ₃ , Sb ₂ O ₃ , others < 0.05% |
| | | | | | |
| | | | | | |



MODERN
Dental Laboratory USA

Heikas Martirosian

Heikas Martirosian, CDT
Technical Director

