

Standard Protocol

The following are our standard protocols that we follow in the event of a question. By giving the lab pre-determined instructions, either by standard protocol or written instructions on the Rx, we can produce quality restorations within the scheduled time frame.

The laboratory reserves the right to use alternative protocol methods, if warranted.


Issue	Standard Protocol
If no occlusal clearance for prescribed restoration	Lab reduces the opposing tooth, marks the opposing model with the reduction.
Occlusal Staining	Light
Metal Design	Porcelain Buccal, Lingual Band
Butt Margin	No
Pontic Design	Modified Ridge Lap
Embrasure	Natural
Request Date	Twelve working days – Max Please use Calendar

Should the laboratory exercise one of the protocols, the case will arrive with a "Doctor's Note" to indicate the procedure the lab followed.

EXAMPLE OF "DOCTOR'S NOTE"



Modern Dental Laboratory USA
877.711.8778
info@moderndentalUSA.com



Creating the finest restorations in the world.

Technical Note to Doctor!

Please advise Doctor of this change prior to seating. Thank you.
See note on reverse.

Case _____

For consultation, please contact our Technical Department.