

# Certificate of Authenticity

The dental restorations used for your case were constructed of FDA approved 510k materials and hand crafted in our North American Digital Processing Center. This is your assurance of quality and safety in materials and construction.

Item	Product Name	Manufacturer	Country of Origin	FDA 510k	Description
Zirconium oxide core	DD Bio ZX <sup>2</sup>	Dental Direkt	Germany	K093748	ZrO <sub>2</sub> , HfO <sub>2</sub> , Y <sub>2</sub> O <sub>3</sub> , Al <sub>2</sub> O <sub>3</sub> , other oxides <0.25%
Porcelain	Vita VM 9 Effect Bonder Fluid	Vita Zahnfabrik	United States	K052710	Feldspathic porcelain (no breakdown given)



*Heikas Martirosian*

Heikas Martirosian, CDT  
Director of Technical Training