

The dental restorations used for your case were constructed of FDA approved 510k materials and hand crafted in our North American Digital Processing Center. This is your assurance of quality and safety in materials and construction.

Item	Product Name	Manufacturer	Country of Origin	FDA 510k	Elements & Compounds/Desc
Esthetic Core	IPS Empress Esthetic	Ivoclar	Germany	K982616	>98% SiO <sub>2</sub> , BaO, Al <sub>2</sub> O <sub>3</sub> , CaO, CeO <sub>2</sub> , NaO <sub>2</sub> , K <sub>2</sub> O, B <sub>2</sub> O <sub>3</sub> , <2% TiO <sub>2</sub> and pigments
Veneer Porcelain	IPS Empress Esthetic Veneer	Ivoclar	Germany	K913372	SiO <sub>2</sub> >55%, Other: Al <sub>2</sub> O <sub>3</sub> , K <sub>2</sub> O, Na <sub>2</sub> O, CaO, B <sub>2</sub> O <sub>3</sub> , TiO <sub>2</sub> , CeO <sub>2</sub> , and pigments



*Heikas Martirosian*

Heikas Martirosian, CDT  
Director of Technical Training