

# Certificate of Authenticity

The dental restorations used for your case were constructed of FDA approved 510k materials and hand crafted in one of our global processing centers. This is your assurance of quality and safety in materials and construction.

| Item               | Product Name             | Manufacturer | Country of Origin | FDA 510k | Description   |
|--------------------|--------------------------|--------------|-------------------|----------|---|
| Alloy, Noble Metal | Callisto CPG             | Ivoclar      | Liechtenstein     | exempt   | Co 39.9%, Pd 24.6%, Cr 21.3%, Mo 11.6%, W 1.2%<br><1% Au, Fe, Ta, B, Al   |
| Porcelain          | GC Initial MC, AL and LF | GC America   | United States     | K040567  | Ceramic Powder made from SiO <sub>2</sub> , Al <sub>2</sub> O <sub>3</sub> , K <sub>2</sub> O, Na <sub>2</sub> O, CaO, BaO, B <sub>2</sub> O <sub>3</sub> , <1% LiO <sub>2</sub> , ZrO <sub>2</sub> , other oxides (CeO <sub>2</sub> , CeF <sub>3</sub> ) |



*Heikas Martirosian*

Heikas Martirosian, CDT  
Director of Technical Training