


## FULL ACRYLIC DENTURE

# Certificate of Authenticity

The dental restorations used for your case were constructed of FDA approved 510k materials and hand crafted in one of our global processing centers. This is your assurance of quality and safety in materials and construction.

Item	Product Name	Manufacturer	Country of Origin	FDA 510k	Description
Standard Acrylic	Vertex™ Regular	Vertex Dental B.V.	Netherlands	K102654	Polymethylmethacrylate and Benzoyl Peroxide; Methyl Methacrylate of ethylene glycol dimethacrylate, N, N-dimethyl-4-toluidine, Tinuvin P
Premium Acrylic	SR Ivoclar High-Impact Acrylic	Ivoclar	Germany	K896130	Polymethylmethacrylate and Benzoyl Peroxide; Methyl Methacrylate, Ethylene glycol dimethacrylate
Standard Acrylic Teeth	Vita MFT®	Vident	Germany	K914476	Highly cross-linked acrylic resin teeth based on polymethylmethacrylate
Premium Denture Teeth	SR-Vivadent & SR-Orthotype	Ivoclar	Germany	K864707	Synthetic polymer teeth based on polymethylmethacrylate



**MODERN**  
Dental Laboratory USA

*Heikas Martirosian*

Heikas Martirosian, CDT  
Director of Technical Training